



Neutron and X-ray computed tomography of a natural uranium tristructural isotropic (TRISO) fuel compact

June 2025

Changing the World's Energy Future

William C Chuirazzi, Swapnil Kishor Morankar, Yuxuan Zhang



DISCLAIMER

This information was prepared as an account of work sponsored by an agency of the U.S. Government. Neither the U.S. Government nor any agency thereof, nor any of their employees, makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness, of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. References herein to any specific commercial product, process, or service by trade name, trade mark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the U.S. Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the U.S. Government or any agency thereof.

Neutron and X-ray computed tomography of a natural uranium tristructural isotropic (TRISO) fuel compact

William C Chuirazzi, Swapnil Kishor Morankar, Yuxuan Zhang

June 2025

**Idaho National Laboratory
Idaho Falls, Idaho 83415**

<http://www.inl.gov>

**Prepared for the
U.S. Department of Energy
Under DOE Idaho Operations Office
Contract DE-AC07-05ID14517**

March 28th, 2025

**William Chuirazzi¹, Yuxuan Zhang²,
and Swapnil Morankar¹**

¹Idaho National Laboratory

²Oak Ridge National Laboratory

Neutron and X-ray computed tomography of a natural uranium tristructural isotropic (TRISO) fuel compact

Battelle Energy Alliance manages INL for the
U.S. Department of Energy's Office of Nuclear Energy






Idaho National Laboratory

Introduction

- Idaho National Laboratory's mission is to **discover, demonstrate** and **secure** innovative **nuclear energy solutions**, other clean energy options and critical infrastructure
- Nondestructive examination techniques support this mission by:
 - Providing information for **modeling** fuel performance
 - Comparing **pre-** and **post-irradiation** effects
 - Informing destructive analysis
 - Saves **time** and **money**
- Quantify material performance in the harsh environment of a reactor to ensure safe commercial operations once deployed
- Ultimately, we aim to **speed up** the **design, development,** and **deployment** of nuclear energy technologies from conceptualization to commercialization



MFC Nuclear Research, Development & Demonstration Capabilities *(with other connected INL capabilities)*



Fabrication

- Experimental Fuel Facilities
- Fuels & Applied Science Building
- Fuel Manufacturing Facility
- Zero Power Physics Reactor
- Analytical Laboratory
- Advanced Fuels Facility

Fresh Fuel Characterization

- Fuels & Applied Science Building
- Analytical Laboratory
- Experimental Fuel Facilities

Irradiation

- Transient Reactor Experiment & Testing (TREAT)
- Neutron Radiography Reactor (NRAD)
- Advanced Test Reactor (ATR)
- Offsite Reactors

Post-Irradiation Examination & Characterization

- Hot Fuel Examination Facility
- Irradiated Materials Characterization Lab
- Fuel Conditioning Facility
- Analytical Laboratory
- Fuels & Applied Science Building
- Electron Microscopy Lab
- Neutron Radiography Reactor (NRAD)

Space Nuclear Power and Isotope Technologies


- Space & Security Power Systems Facility
- Engineering Development Lab
- Idaho Nuclear Technology & Engineering Center

Advanced Reactor Demonstration Test Beds


- TREAT micro-Reactor Experiment Cell
- Laboratory for Operations and Testing in the US
- Demonstration of Microreactor Experiments

TRISO Fuel and Compacts

- TRISO- TRi-structural ISOtropic particle fuel
- PIE of individual particles investigates changes in individual particles layers
 - Thickness changes, buffer delamination, etc.
 - Kernel porosity and fission product migration
- PIE XCT of intact fuel compacts and loose particles yields different information. Fuel compact examination provides:
 - Kernel morphology
 - Variations morphology as a function of spatial position
 - Quantification of the buffer fracture frequency
 - Preservers compact for additional exams
 - Impossible following traditional examinations



ZEISS Xradia 620 Versa at IMCL



- Best spatial resolution to date: **270 nm/voxel**
- Hottest sample to date: **120 R/hr on contact**
- Instrument is maintained as **radiologically clean**


X-ray Source
Micro-focus
Tungsten Target
30 kVp - 160 kVp

Optics & Detectors
Flat Panel Detector
30.7 cm x 19.4 cm Max FOV
100 - 6 μm pixel size




Thermomechanical Load Stage
5000 N Max Load
-20 °C - +160 °C Temp Range
Compression & Tension


Optics & Detectors
4 Objective Lenses
CCD Camera
0.6 mm x 0.6 mm Min FOV
55 - 0.2 μm pixel size






AGR TRISO Fuel




2D slices of a 3D reconstructed volume showing a) a low-energy (40 keV) scan on an irradiated AGR-2 TRISO particle, b) the corresponding high-energy (110 keV) scan of the same particle, c) the images fused together, and d) a 3D rendering of the particle using information from both datasets.



TRISO Compact Examination





XCT Renderings





Quantitative Results

- Quantify kernel information
 - Equivalent diameter
 - Sphericity
 - Kernel position
 - Packing fraction
- Useful to compare pre- and post-irradiation to examine changes
 - Kernel swelling
 - Aspherical particles





Need for Neutron CT

- Can we see defects in the matrix of the compacts nondestructively?
- **Informs structural integrity** of the compact before and after irradiation
- If features of interest (i.e. cracks) are found, their location can be used to determine potential causes
- X-ray attenuation introduces too much noise for matrix examination
 - Advanced reconstruction and artifact reduction techniques cannot sufficiently resolve features in the matrix
- Neutrons may offer a solution
- Used an unirradiated, natural Uranium compact as a proof-of-principal




MARS Beamline

- Multimodal Advanced Radiography Station at HFIR's CG-1D Beamline
- Polychromatic beam of cold neutrons (peak wavelength 2.6 Å)




MARS Imaging Conditions

- Absolute neutron flux was $\sim 5 \times 10^6$ n/cm²/s
- Pixel size is ~ 16 μ m/pixel
 - Spatial resolution is ~ 48 μ m
- L/D ~ 600
- sCMOS camera
- 1091 Radiographs captured
- SIRT Reconstruction algorithm
- Ring artifact removal but otherwise no image processing
 - bm3d ring removal






nCT Results

- Difficult to get information about individual particles
- **Crack** within matrix can **clearly** be **seen** across multiple slices
- Can provide useful method for matrix inspection




nCT Results






nCT Renderings



11







Crack Identification



Future Capabilities – nCT for Multiple FOV's

- Different applications require different FOV's
- A multi-FOV imaging system would cover all applications possible in the NRS beam
- Targeting $\sim 20 \mu\text{m}$ effective spatial resolution ($7 \mu\text{m}$ pixel)



batteries, biofuels)



– Model V&V using realistic geometry

Conclusions

- Neutron CT successfully identified cracking in the graphite matrix
- X-ray and neutron CT perform complementary examinations of fuel compacts
 - **XCT**: Kernel information
 - **nCT**: Matrix information

Potential Future Work:

- nCT of irradiated compact
- XCT/nCT data registration
- Isotopic fission product mapping (nGI/ToF)



Acknowledgements

- Thanks to Phil Winston, Dave Laug, BJ Camphouse, and John Stempien for programmatic support and resources in providing samples and coordinating sample transfer between INL and ORNL.
- Thanks to the HFIR support staff for facilitating the experiment.
- This work was supported by Battelle Energy Alliance, LLC, under Contract No. DE-AC07–05ID14517 with the U.S. Department of Energy.
- This research used resources at the High Flux Isotope Reactor a DOE Office of Science User Facility operated by the Oak Ridge National Laboratory. The beam time was allocated to MARS (CG-1D) on proposal number IPTS-34736.



Idaho National Laboratory

Battelle Energy Alliance manages INL for the U.S. Department of Energy's Office of Nuclear Energy. INL is the nation's center for nuclear energy research and development, and also performs research in each of DOE's strategic goal areas: energy, national security, science and the environment.