

# **2020 Winter Conference on Plasma Spectrochemistry - Luiza Albuquerque - Poster**

Luiza Gimenes Rodrigues Albuquerque,  
Beau J Barker, Leah N Squires

January 2020



The INL is a U.S. Department of Energy National Laboratory  
operated by Battelle Energy Alliance

## **2020 Winter Conference on Plasma Spectrochemistry - Luiza Albuquerque - Poster**

**Luiza Gimenes Rodrigues Albuquerque, Beau J Barker, Leah N Squires**

**January 2020**

**Idaho National Laboratory  
Idaho Falls, Idaho 83415**

**<http://www.inl.gov>**

**Prepared for the  
U.S. Department of Energy**

**Under DOE Idaho Operations Office  
Contract DE-AC07-05ID14517**



# ASSESSMENT OF IMPURITIES IN Np METAL USING DOUBLE-FOCUSING SF-ICP-MS

Luiza G. R. Albuquerque<sup>1\*</sup>; Beau J. Barker<sup>1</sup>; Leah N. Squires<sup>2</sup>

\*luiza.albuquerque@inl.gov

<sup>1</sup>Analytical Laboratory – Materials and Fuels Complex, Idaho National Laboratory – P.O. Box 1625, Idaho Falls – ID 83415

<sup>2</sup>Advanced Fuel Manufacturing and Development – Materials and Fuels Complex, Idaho National Laboratory – P.O. Box 1625, Idaho Falls – ID 83415

**Production of 99.999% pure Np metal from NpO<sub>2</sub> and the determination of its impurities are the main goals of this research**

- High-purity Np is required to accurately measure its fundamental physical properties, in order to support transmutation and nuclear fuels development research.
- Very little Np is available in its metal form, so a process of producing Np metal from NpO<sub>2</sub> was recently developed at INL<sup>1</sup>.


[1] Squires, L. and Lessing, P. 2016. *Journal of Nuclear Materials*, 471, pp. 65-68  
[2] Richter, S. et al., 2013. *Journal of Analytical Atomic Spectrometry*, 28, 10, pp. 1540-1543  
[3] Becker, J. S. and Dietze, H. 2003. *International Journal of Mass Spectrometry*, 228, 2-3, pp. 127-150

**Purity of the material (%) = 100% Np – Impurities<sup>2,3</sup>**  
**Maximum certifiable purity = 100% Np -  $\sum$  LOQ**

Analytical challenges:

***Identify impurities and achieve low detection limits***

**NpO<sub>2</sub> reduction**



**A neptunium metal button is formed at the bottom of the reaction mixture due to its higher density**

**Instrumentation**

## Instrumental parameters

RF Power (W)	1300
Ar cooling gas ( $L \text{ min}^{-1}$ )	14
Ar auxiliary gas ( $L \text{ min}^{-1}$ )	1
Nebulizer gas pressure (psi)	32
Sample uptake rate ( $\mu L \text{ min}^{-1}$ )	190
Dwell time per peak (ms)	5
Number of sweeps	300
Resolution modes ( $m/\Delta m$ )	300 and 4000
Acquisition modes	300 RP Deflector jump 4000 RP Deflector scan




Nu Instruments Attom ES HR-ICP-MS  
Attached to a radiological hood

**Sample dissolution**




Np metal button  
(approx. 0.6 g)

Sequential washes  
removal of surface  
contaminants




Dissolution  
25 mL 6M HCl  
+ 150  $\mu L$  24M HF



90 minutes

Matrix conversion to HNO<sub>3</sub>  
10 mL 16M HNO<sub>3</sub>



80 minutes


**Final solution  
50 mL 8M HNO<sub>3</sub>**



Both washes and the dissolved sample were first evaluated by its qualitative mass scans. Elements were later quantified using Double Focusing SF-ICP-MS

**Impurities identification**

- The water wash removed mostly analytes present in the Np metal production process, indicating that they were present as surface contamination.
- The HCl wash had a similar peak pattern to the Np metal sample, indicating that the process not only removed surface contamination, but also reacted with the sample itself.
- Tantalum, and some actinides are the candidates for major impurities on the Np metal. Minor impurities could include Al, Au, Ca, Cr, Fe, I, Nb, Ni, Pb, Sn, Tl, Y and Zn. These contaminants can have several different origins, such as the materials used for the casting process, the use of contaminated NpO<sub>2</sub> as the starting material or environmental contamination.





# ASSESSMENT OF IMPURITIES IN Np METAL USING DOUBLE-FOCUSING SF-ICP-MS


*Luiza G. R. Albuquerque<sup>1\*</sup>; Beau J. Barker<sup>1</sup>; Leah N. Squires<sup>2</sup>*

*\*luiza.albuquerque@inl.gov*

<sup>1</sup>Analytical Laboratory – Materials and Fuels Complex, Idaho National Laboratory – P.O. Box 1625, Idaho Falls – ID 83415


<sup>2</sup>Advanced Fuel Manufacturing and Development – Materials and Fuels Complex, Idaho National Laboratory – P.O. Box 1625, Idaho Falls – ID 83415

Impurities evaluation and HR-IP-MS analyses



Schematic of the Attom ES: High Resolution ICP-MS (Nu Instruments)

- Double-focusing instrument – ESA focuses the energy of the ions coming from the entrance slit, and magnet focuses the ions related to their m/z onto the exit slit
- Deflectors voltage can be changed for each “parked mass” of the magnet – faster analysis over a wide mass range
- Adjustable entrance and exit slits are used to set the Resolving Power (RP), varying from 300RP to 10000RP
- A narrower slit increases the RP, but decreases sensitivity due to attenuation of the ion beam



Mass spectrum of a 0.2 ng g<sup>-1</sup> Cr solution


## Example: <sup>52</sup>Cr<sup>+</sup> evaluation in Np samples

Resolving Power	300	4000
5% HNO <sub>3</sub> signal intensity (cps)	369000	1820
Calibration slope	729000	121000
Calibration r <sup>2</sup>	0.9957	0.9964
LOQ (μg g <sup>-1</sup> )	40	0.3
Measurement RSD (%)	5 – 20	1 – 15
<sup>52</sup> Cr <sup>+</sup> concentration (μg g <sup>-1</sup> )	<40	34.3 ± 1

$$LOQ = 15 * \frac{SD \text{ of } 10 \text{ measurements of the blank}}{\text{Slope}} * \text{sample dilution factor}$$

- Background signal reduction of over 2 orders of magnitude  
*Interference from <sup>40</sup>Ar<sup>12</sup>C<sup>+</sup> resolved*
- Sensitivity reduced by a factor of 5  
*As expected, due to beam attenuation*
- No significant change in determination coefficient
- Quantification limit increased by 2 orders of magnitude  
*Eliminating interference compensates for sensitivity loss*
- Slight improvement
- <sup>52</sup>Cr could be quantified using a higher RP

- Tantalum and Fe were the main impurities, being introduced in the sample from Ta stirrer and the stainless steel crucible used in the casting process.
- Uranium, Th and Pu could have originated from an impure starting NpO<sub>2</sub> material.
- Calcium was also one of the contributors, originating from the Ca/CaCl<sub>2</sub> used in the reduction process



Isotope/emission line used						
Element/Isotope	RP	for total conc. calculation	Interferences	Conc. (μg/g)	RSD (%)	LOQ (μg/g)
Al	4000	<sup>27</sup> Al	-	28.9	±3	3
Cr	4000	<sup>52</sup> Cr	<sup>36</sup> Ar <sup>16</sup> O <sup>+</sup> , <sup>38</sup> Ar <sup>14</sup> N <sup>+</sup> , <sup>40</sup> Ar <sup>12</sup> C <sup>+</sup>	41	±1	0.4
Fe	4000	<sup>56</sup> Fe	<sup>40</sup> Ar <sup>16</sup> O <sup>+</sup>	396	±2	50
I	4000	<sup>127</sup> I	<sup>238</sup> U <sup>16</sup> O <sup>++</sup>	<4	N/A	4
Nb	300	<sup>93</sup> Nb	-	1.55	±4	0.003
Ni	300	<sup>60</sup> Ni	-	19.6	±7	0.8
Pb	300	<sup>208</sup> Pb	-	<0.1	N/A	0.1
Sn	4000	<sup>118</sup> Sn	-	<0.2	N/A	0.2
Ta	300	<sup>181</sup> Ta	-	1920	±2	0.1
Th	300	<sup>232</sup> Th	-	103	±0.1	0.006
Tl	300	<sup>203</sup> Tl	-	<0.001	N/A	0.001
V	4000	<sup>51</sup> V	<sup>14</sup> N <sup>37</sup> Cl <sup>+</sup>	1.68	±6	0.02
Y	300	<sup>89</sup> Y	-	10.8	±1	0.005
<sup>233</sup> Pa	300	<sup>233</sup> Pa	-	4.56	±7	0.7
<sup>234</sup> U	300	<sup>234</sup> U	-	9.53	±1	0.001
<sup>235</sup> U	300	<sup>235</sup> U	-	2.74	±1	0.001
<sup>236</sup> U	300	<sup>236</sup> U	-	4.25	±5	0.5
<sup>238</sup> U	300	<sup>238</sup> U	-	115	±1	0.01
<sup>239</sup> Pu	300	<sup>239</sup> Pu	-	35.6	±1	0.002
<sup>240</sup> Pu	300	<sup>240</sup> Pu	-	6.05	±1	0.002
<sup>241</sup> Am	300	<sup>241</sup> Am	-	1.16	±7	0.002
<sup>242</sup> Cm	300	<sup>242</sup> Cm	-	2.55	±2	0.001
<sup>243</sup> Am	300	<sup>243</sup> Am	-	0.243	±5	0.001
<sup>244</sup> Cm	300	<sup>244</sup> Cm	-	0.0837	±4	0.001
Ca <sup>+</sup>		393.366 nm		97	±1	20

\* Calcium was analyzed by ICP OES

Total impurities (μg g<sup>-1</sup>)  
Total Np purity (%)  
Maximum certifiable purity (% as 100% Np - Σ LOQ)

2800 ± 2%  
99.720  
99.992

- ❖ Np metal that was produced did not achieve the desired purity (99.999%)
- ❖ It was possible to identify the main sources of impurities, with great precision and accuracy
  - The casting process will be modified accordingly to the results obtained for the main contaminants.
- ❖ Using the current detection limits, we could certify a material with up to 99.992% purity
  - Detection limits can be improved - using of intermediate RP (2500).
  - Could also be improved by removing the Np from the solution allowing a more concentrated solution to be used and removing potential spectral interferences from the ICP OES measurements.

Acknowledgements

This work was supported by the U.S. Department of Energy, Office of Nuclear Energy under DOE Idaho Operations Office Contract DE-AC07- 051D14517 as part of a Nuclear Science User Facilities project



Dr. Vivian M. de O. Carioni  
HR-ICP-MS Instrument Scientist

Scott R. Wilde  
Np Metal Reduction

Cortney M. Pincock  
Q-ICP-MS Instrument Scientist